

LID FOR MALE URINAL BULBOUS

Eco-Friendly Medical Pulp



Made from
BIODEGRADABLE
Double Lined Kraft

Designed for Single Use to
Reduce Risk of Infection

The Lid for Male Urinal Bulbous is made from eco-friendly biodegradable kraft material, intended for use as a cap fitted over the opening of male urinal bulbous to prevent spillage and reduce odour.



Manufactured as single-use disposable medical pulp to meet the requirement to hold water and/or human waste for a minimum period of four hours at an initial temperature of 35°C. Suitable for disposal into macerator which complies with PAS29:1999 "Disposable Pulp For Use in Healthcare".

Suitable for



Hospitals



Nursing Homes



Personal Use




Medical Clinics




Disability Homes


User Guide




1 Remove the lid from Male Urinal Bulbous and before use.



2 Place the Male Urinal Bulbous between the legs for urination.



3a Dispose the Male Urinal Bulbous together with the lid directly into macerator or empty the content from the Male Urinal Bulbous before disposing into appropriate waste bins.



Stackable



Water Retention



Macerator Friendly



Eco-Friendly



Recycled

PRODUCT DESCRIPTION

Adaptor for Urinal Bulbous

PRODUCT SIZE

85 x 70 x 15 mm

PACKAGING

1000 pcs per carton

CARTON DIMENSION

510 x 260 x 360 mm

ORDER CODE

100APPLMU1000



Manufactured in accordance with PAS29:1999
Disposable Pulp Products For Use in Healthcare

Rev2(0225)EN



APEX PULP PTE LTD

www.apexpulp.com

sales@apexpulp.com

© 2025 APEX PULP PTE LTD. All rights reserved.