

INJECTION TRAY 200ML Eco-Friendly Medical Pulp

Made from BIODEGRADABLE Sugarcane Bagasse


Designed for Single Use to Reduce Risk of Infection

© 2025 APEX PULP PTE LTD. All rights reserved.






User Guide



Place surgical items on the Injection Tray.



After use, dispose the Injection Tray into macerator or appropriate waste bins

Injection Tray 200ml Eco-Friendly Medical Pulp

The Food Grade Injection Tray with 200ml rated capacity is made from eco-friendly biodegradable sugarcane bagasse, intended for use as multi-purpose tray to carry syringes and other medical items. It is designed with smooth inner finishing and meshed outer surface for extra grip.

Manufactured as single-use disposable medical pulp to meet the requirement to hold water and/or human waste for a minimum period of four hours at an initial temperature of 35°C. Suitable for disposal into macerator which complies with PAS29:1999 "Disposable Pulp For Use in Healthcare".

Features

- Ideal for surgical procedures
- Able to retain water for at least 4 hours
- Non-sterile, single use tray
- Meshed outer surface for better grip
- Made from sugarcane bagasse
- Biodegradable and compostable
- Can be disposed into macerator
- Manufactured in accordance with:
 - EU Regulation No. 1935/2004 Materials intended to come into contact with food
 - US FDA 21 CFR 176.170 Components of paper and paperboard in contact with aqueous and fatty foods.

Caution

• Do not pour alcohol based solutions into the Injection Tray.





Water Retention



Macerator Friendly



Eco-Friendlu

PRODUCT DESCRIPTION	PRODUCT SIZE	PACKAGING	CARTON DIMENSION	ORDER CODE
Injection Tray 200ml	200 x 75 x 30 mm	1000 pcs per carton	400 x 300 x 230 mm	101APPTIJ0200
)E N





Manufactured in accordance with PAS29:1999 Disposable Pulp Products For Use in Healthcare



APEX PULP PTE LTD www.apexpulp.com

sales@apexpulp.com © 2025 APEX PULP PTE LTD. All rights reserved.